

Hilfsmittel

Single Use Bacterial/Viral Filter



Product Feature

- Machine End Device
- Bacterial / Viral
- 99.999% Efficiency
- Electrostatic Media
- 1/4" ISO connectors
- Available in blue or green ring

Cat. No.	Description	Pack Per Box
CM-AO-1001	Suction/Inflution Filter	50
	Clinical Indication	Suction/insufflation filter
	Filtration Efficiency	99.999% Efficiency
	Recommended Usage	Up to 48 hours
	Casing Material	Polypropylene

Comcorde Medical (M) Sdn. Bhd.

No. 1, Jalan 11/32B,
TSI Business Industrial Park,
6½ Mile Off Jalan Kepong,
52000 Kuala Lumpur
Malaysia

Tel No 603-6250 9657
Fax No 603-6250 5367
Website www.comcorde.com.my
Email sales@comcorde.com.my

