


Hilfsmittel®

Pediatric Bacterial / Viral Filter (HME)



Product Feature

- Anesthesia/Ventilation
- Clear Casing
- Electrostatic Media
- Bacterial/viral
- 22M / 15F-15M ISO connectors
- Transparent housing
- WITH gas sampling port only
- Individual packed, available in sterile or non-sterile

Cat. No.	Description	Pack Per Box																										
CM-AO-9066	Pediatric HME filter WITH gas sampling port	50																										
	<table border="1"> <tr> <td>Clinical Indication</td> <td>Short Term Anesthesia</td> </tr> <tr> <td>Tidal Volume Range</td> <td>≥ 75mL</td> </tr> <tr> <td>Penetration</td> <td>At 15 L.min⁻¹ (mean [range]) Unused 0.52[0.49to0.54] % After 3h simulated use 0.55[0.53to0.58] %</td> </tr> <tr> <td>Pressure Drop</td> <td>At 15 L.min⁻¹ Unused 70Pa Increase after 3h simulated use 1Pa</td> </tr> <tr> <td>Moisture Output</td> <td>0.25 L tidal volume 3 lgm⁻³</td> </tr> <tr> <td>Weight</td> <td>17.5g</td> </tr> <tr> <td>Dimensions (Diameter x L)</td> <td>52 x 66 mm</td> </tr> <tr> <td>Dead Space</td> <td>33 mL</td> </tr> <tr> <td>Connector</td> <td>Compliance with BS EN 1281-1 Patient End 15F /22M Machine End 15M</td> </tr> <tr> <td>Filtration Efficiency</td> <td>99.999% Efficiency</td> </tr> <tr> <td>Casing Material</td> <td>Styrolux</td> </tr> <tr> <td>HME Material</td> <td>Polyurethane Foam treated with calcium Chloride</td> </tr> <tr> <td>Recommended Usage</td> <td>Up to 48 hours</td> </tr> </table>	Clinical Indication	Short Term Anesthesia	Tidal Volume Range	≥ 75mL	Penetration	At 15 L.min ⁻¹ (mean [range]) Unused 0.52[0.49to0.54] % After 3h simulated use 0.55[0.53to0.58] %	Pressure Drop	At 15 L.min ⁻¹ Unused 70Pa Increase after 3h simulated use 1Pa	Moisture Output	0.25 L tidal volume 3 lgm ⁻³	Weight	17.5g	Dimensions (Diameter x L)	52 x 66 mm	Dead Space	33 mL	Connector	Compliance with BS EN 1281-1 Patient End 15F /22M Machine End 15M	Filtration Efficiency	99.999% Efficiency	Casing Material	Styrolux	HME Material	Polyurethane Foam treated with calcium Chloride	Recommended Usage	Up to 48 hours	
Clinical Indication	Short Term Anesthesia																											
Tidal Volume Range	≥ 75mL																											
Penetration	At 15 L.min ⁻¹ (mean [range]) Unused 0.52[0.49to0.54] % After 3h simulated use 0.55[0.53to0.58] %																											
Pressure Drop	At 15 L.min ⁻¹ Unused 70Pa Increase after 3h simulated use 1Pa																											
Moisture Output	0.25 L tidal volume 3 lgm ⁻³																											
Weight	17.5g																											
Dimensions (Diameter x L)	52 x 66 mm																											
Dead Space	33 mL																											
Connector	Compliance with BS EN 1281-1 Patient End 15F /22M Machine End 15M																											
Filtration Efficiency	99.999% Efficiency																											
Casing Material	Styrolux																											
HME Material	Polyurethane Foam treated with calcium Chloride																											
Recommended Usage	Up to 48 hours																											

CoMcorde Medical (M) Sdn. Bhd.

No. 1, Jalan 11/32B,
TSI Business Industrial Park,
6½ Mile Off Jalan Kepong,
52000 Kuala Lumpur
Malaysia



Tel No 603-62433245

Fax No 603-62433246

Website

www.comcorde.com.my

Email

sales@comcorde.com.my

